



Predator attack Canadian mountainous terrain, likely bear or cougar

Pregnant brood mare, Alberta, Canada

Multi Radiance ACTIVet 25W Super Pulse Laser

No pain meds after week one (1).

Antibiotics as precaution.

Fly mask and dressing used daily (head bandage for eye protection and to securely hold gauze in place)

Cleaned wound gently daily with water only.

Treated 5 days/week for 2 weeks (twice daily)

Head wound was completely closed, scratches healed and eyesight saved with normal function with peripheral vision intact

Days 1–30 treatment progression

- 1st day : Blue followed by 1000-3000Hz 5 min
- 2nd treatment: PHT 50Hz over wound + 1-250 Hz 30 seconds... drainage site kept open
- Swollen eye lid - hematoma likely without drainage 1000 -3000 Hz plus Blue on eyelid (only) for pain relief and swelling reduction
- Week 3 - 4 one treatment/day only
- Week 4 - 3 times (once treatment/day only)
- Returned to herd after 3rd week, fly mask only



© Lisa Cawte-Baker



© Lisa Cawte-Baker